


## The first stent CE Marked for endoscopic ultrasound guided transluminal therapy


Cautery-enabled access and delivery catheter with the preloaded therapeutic AXIOS Stent for an exchange-free procedure.

- Proprietary one-step combined diathermic ring and cut-wire provides easy access into target tissue
- Perpendicular flanges secure tissue layers and help to prevent migration
- Stent creates anastomotic conduit between two lumens
- Large diameter fully covered lumen apposing stent enables drainage and passage of the endoscope for additional therapeutic procedures





Indicated for use to facilitate transgastric or transduodenal endoscopic drainage of a pancreatic pseudocyst, walled off necrosis (>/= 70% fluid content) or the biliary tract.\*




## Ordering Information

## Hot AXIOS™

Stent and Electrocautery Enhanced Delivery System

Order Number	Lumen Diameter (mm)	Saddle Length (mm)	Flange Diameter (mm)	Catheter Outer Diameter (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)
M005 <b>5352</b> 0	6	8	14	9	138	144
M005 <b>5353</b> 0	8	8	17	9	138	144
M005 <b>5354</b> 0	10	10	21	10.8	138	144
M005 <b>5355</b> 0	15	10	24	10.8	138	144
M005 <b>5356</b> 0	20	10	29	10.8	138	144







WARNING: The safety and effectiveness of the 20mm Stent for drainage of the Gallbladder and Bile Duct has not been established.

Contact your local Boston Scientific representative for more information.

